

Product Parameter

Front wheel	8" wheel
Rear Wheel	12.5" wheel
Net Weight	26.8kg
Brake	Electromagnetic Braking system
Batteries	Lithium Battery 24V 11Ah
Motors (pair)	Brushless Motor 180W
Charger	High Power 24V 2A
Controller	Reliable Brushless Controller
Maximum Range	15-18 km
Suspension System	Yes
Ground Clearance	8 cm
Max Speed	8 km/h
Charging Hour	5-6 Hours
Max Slope	12 degrees
Weight Capacity	150 KG
Dimension (Folded)	36*60*78cm
Dimension (Folding)	108*60*99cm

1. Range and speed vary with user weight terrain type, battery charge battery condition and tyre pressure.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
3. Battery weight may vary based on manufacturer.

TO BE A YATTL AGENT

2021

Requirements

1. At least 1 lawful and registered local shop or sale place managed by yourself.
2. Legally registered company.
3. Have some knowledge and strong interest in electric wheelchairs.

Our Support

1. Agency area protection.
2. The best quality and competitive price.
3. Prior rights to the new products with new opportunity.
4. No minimum ordering quantity.
5. Soonest response on any of your questions.
6. Free parts supply support policy.

JOIN US



YATTL INDUSTRY CO., LTD.



Tel: +86 755 8170 7199

Web: www.yattll.com

Email: info@yattll.com





Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China.

Office: RM.1403, BLK.B of Duocai Office Complex, Gu-anlan High-tech Zone, Longhua DIST., Shenzhen, China.

New upgrade
**Reinforce
frame**

Model: **YE100**



www.yattll.com |    





YE100 Lorem Ipsum

YE100 is with an upgraded frame base on the original models, which makes the chair stronger and more stayble. The maximum loading capacity increase to 150kg. The new designed seat cushion is more fashionable and comfortable.

Features

1. Redesigned & Reinforced frame with suspension on back frame.
2. Brushless motor 24V 180W each.
3. New designed cushion, comfortable and removable.
4. SAMSUNG lithium battery 24V 11Ah.
5. 8" PU Solid front wheel, 12.5" back wheel.
6. Quick release joystick.
7. Foldable footrest to make passengers be closer to the seat.

Colors

Cushion color




Frame color



Folding steps




Three steps quick folding , Only takes 1s.




Suspension system



Foldable legrest



Reinforced frame



Intelligent joystick

 • YATTL INDUSTRY CO., LTD.

Tel: +86 755 8170 7199

Email: info@yattll.com