

Economical Model **YE145D**



New Design

Product Parameters

Model No. YE145D	Recline NO
Seat width 45.5 cm	Seat depth 50 cm
Front wheel 8" Suspension wheel	Rear wheel 12.5" Solid, Suspension wheel
Battery Lithium 11Ah	Motor Brushed 250W each
Charger 24V 2A	Controller Brushed controller
Brakes Electromagnetic braking	Drive distance 12-15 km
Max speed 6 km/h	Net weight 25 kg
Weight capacity 135 kg	Expanded Size 107*61*100cm

EASY TRAVELLING



Main Features:

- Ergonomic- designed seatback cushion for more comforts
- Aircraft aluminium alloy frame
- Minimaft but not simple
- Unique wheels with perforations
- Intelligent controls 360° rotation




**Powerful brushed motors
250W*2**



**Intelligent controller
Easy control**



**Suspension wheel
Rubber material**



One step fold