

Reinforced Frame

YE200



Main Features:

- Unique-designed battery assembled on Two sides
- High precision craftsmanship
- Equipped with brushless motors 250W*2
- Ergonomic& highly comfortable seat back cushion
- Loading capacity 150kg
- Driving range 15-18km



Unique Dual-battery Design



Brushless motors 250W*2



High-precision craftsmanship



12.5" PU solid tyre



Product parameters

Model No. YE200	Recline NO
Seat width 49 / 52 / 56 cm	Seat depth 43 cm
Front wheel 8"PU solid wheel	Rear wheel 12.5"PU solid wheel
Battery 6Ah/9Ah *2pcs	Motor Brushless 250W
Charger 24V 2A/3A/5A	Controller Brushless controller
Brakes Electromagnetic braking	Drive distance 15-18 km
Max speed 6 km/h	Net weight 25.5 kg
Weight capacity 150 kg	Expanded Size 60*98*110 cm