

Carbon Fiber Model

YECX1



Carbon Fiber 12.5KG


- Weighing 12.5KG w/o battery.
- Magnesium alloy hub wheel.
- Brushless motor 24V 300W.
- Optional lithium battery 6Ah/10Ah/20Ah.

Product Parameters

Seat width	46 cm
Seat depth	41 cm
Front wheel	7" PU
Rear wheel	12" PU
Battery	6Ah / 10Ah / 20Ah
Motor	Brushless motor 300W
Charger	24V 2A
Controller	Intelligent brushless controller
Brakes	Electromagnetic braking
Drive distance	15-17 km
Weight capacity	100 kg
Max speed	6 km/h
Net weight	12.5 kg
Expanded Size	900*575*870 mm



REBOOT
YOUR
LIFE



Carbon Fiber Frame

Ultra-lightweight carbon fiber material T700,
net weight of 12.5kg only.



Intelligent LCD Controler

User-Friendly, One-Hand Control

