



YE245CR ELECTRIC BACKREST RECLINING

Foldable & Backrest Adjustable

N.W. (W/O BAT)

25.5KG

WEIGHT CAPACITY

135 KG

RECLINING ANGLE

125°



Recline.

Lift up.



Net weight
27.3kg

Seat width
45 cm

Seat depth
45 cm

Reclining Angle
125°

Brake
Electromagnetic brake

Motor
Brushless motor 180W

Battery
Lithium battery 24V 11Ah

Controller
Programmable Controllers

Front wheel size
8 inches

Rear wheel size
12.5 inches

charger
24V/2A

Obstacles
4 cm

weight capacity
135 kg

Maximum speed
6 km/h

Maximum climb
12°

Maximum mileage
15-18 km

Unfolded size
103*62*98 cm

Folded size
49*62*84 cm

TO BE A YATTL AGENT

JOIN US

REQUIREMENTS

1. At least 1 lawful and registered local shop or sale place managed by yourself.
2. Legally registered company.
3. Have some knowledge and strong interest in electric wheelchairs.



OUR SUPPORT

1. Agency area protection.
2. The best quality and competitive price.
3. Prior rights to the new products with new opportunity.
4. No minimum ordering quantity.
5. Soonest response on any of your questions.
6. Free parts supply support policy.


Details display chart



Brushless motor



Angle & Height adjustable



Lift motor & Air cylinder

NEW PRODUCT

- This wheelchair is designed with an awesome Reclining backrest, and comfortable leg rest.
- It has a medical motor and a air cylinder, which makes the reclining and lifting very smooth, and holds 150kg at least.
- Two steps to fold and unfold, making it very simple to use for all passengers.

Features

1. Reclining back rest and Extendable leg rest.
2. Angle and Height adjustable Head rest.
3. Easily foldable and removable.
4. Maximum loading capacity: 135KG.

