

# T1000 CARBON FIBER

YECZ1



## Product Parameters

Model No. YECZ1	Net weight 13 kg
Seat width 45 cm	Seat depth 40 cm
Front wheel 7" inches	Rear wheel 8" inches
Battery Lithium 10.4Ah	Motor Brushless 180W each
Charger 24V 2A	Controller Brushless controller
Max speed 6 km/h	Drive distance 16-20 km
Weight capacity 125 kg	Expanded size 975 * 605 * 945 mm

# JUST 13KG ULTRA-LIGHT FRAME



Net weight  
13 KG  
One-touch folding

Loading capacity  
125kg



Brushless motor  
180W



## Product Details /

01

Smart control, flexible options, with or without display.



02

Plug and charge — easy and convenient

03

Made of T1000 carbon fiber — ultra-light at only 13kg.



04

Wear-resistant Tires, Suitable for All Terrains

Email: [info@yattll.com](mailto:info@yattll.com)

Factory: YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd, Tangxia, Dongguan City, Guangdong, China.

Office: Rm 907, Yichuang International Centre, No. 8 Gaoerfu Blvd, Longhua Dist, Shenzhen, Guangdong, China



YATLL INDUSTRY CO.,LTD

755 - 817 - 07199 [www.yattll.com](http://www.yattll.com)

