

YATTLL YE135

Stylish Purple, Stand Out in Every Journey



Lightweight Strength, Freedom to Go



Drive distance	15-18 km
Seat width	45cm
Seat depth	43 cm
Front wheel	7" solid wheel
Front wheel	12" solid wheel
Weight capacity	135 kg
Max speed	6 km/h

Battery	Lithium battery 11Ah
Motor	Brushless 180W *2
Charger	24V 2A
Controller	Brushless controller
Brakes	Electromagnetic braking
Net weight	22.5 kg
Expanded Size	84*40*60 cm

YE135

4-SPRING SUSPENSION

Smoother Every Journey

