

T1000 CARBON FIBER

YECZ1



Product Parameters



Model No. YECZ1	Net weight 13 kg
Seat width 45 cm	Seat depth 40 cm
Front wheel 7" inches	Rear wheel 8" inches
Battery 11Ah / 15Ah / 20Ah	Motor Brushless 180W each
Charger 24V 2A	Controller Brushless controller
Max speed 6 km/h	Drive distance 16-25 km
Weight capacity 125 kg	Expanded size 975 * 605 * 945 mm

JUST 13KG ULTRA-LIGHT FRAME



Net weight
13 KG
One-touch folding

Loading capacity
125kg



Brushless motor
180W

Product Details /

01

Smart control, flexible options , with or without display.



02

Foldable footrest, freely configurable.

03

Made of T1000 carbon fiber — ultra-light at only 13kg.



04

Plug and charge — easy and convenient.

Email: info@yattll.com
Factory: YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd, Tangxia, Dongguan City, Guangdong, China.
Office: Rm 907, Yichuang International Centre, No. 8 Gaoerfu blvd, Longhua Dist, Shenzhen, Guangdong, China

YATLL INDUSTRY CO.,LTD
755 - 817 - 07199 www.yattll.com

