

ALL CARBON FIBER SCOOTER

“

Carbon Fiber · Only 11 kg

”



11 kg

Net weight

180 W

Brushless motors

125 kg

Loading capacity



INTEGRATED CARBON FIBER TILLER

Lightweight yet strong, with key ignition and intuitive controls placed within easy reach.

LIFT-AND-GO BATTERY

Two removable batteries provide flexible charging and greater freedom for every journey.



Easy to Carry · One-Step Folding



Drive distance	16–20 km	Front wheel	7 inches	Battery	6Ah * 2pcs	Net weight	11 kg
Seat width	38 cm	Rear wheel	7 inches	Motor	Brushless 180w	Weight capacity	125 kg
Seat depth	30 cm	Expanded size	100*45*87 cm	Charger	24V 2A		
Max speed	6 km/h	Fold size	100*45*32 cm	Brakes	Electromagnetic braking		

