

PRODUCT DATA

TECHNICAL SPECIFICATIONS

MODEL	YECZ4 A
SEAT WIDTH	468 MM
SEAT DEPTH	320 MM
FRONT WHEEL SIZE	7 INCH
REAR WHEEL SIZE	7 INCH
BATTERY CAPACITY	6 AH × 2
MOTOR POWER	200W
CHARGER	24V 2A
MAXIMUM LOAD CAPACITY	180KG
MAXIMUM SPEED	6 KM/H
NET WEIGHT	14KG
UNFOLDED (L × W × H)	992 × 468 × 887 MM

Tel: +86 755 81707199

Email: info@yattll.com

Website: www.yattll.com

Factory: YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd,
Tangxia, Dongguan City, Guangdong, China.

Office: Rm 907, Yichuang International Center, No. 8 Gaoerfu Blvd,
Longhua Dist, Shenzhen, Guangdong, China



YATLL INDUSTRY CO.,LTD

PRODUCT BROCHURE

Z4A

CARBON FIBER MOBILITY SCOOTER BUILT FOR EVERY JOURNEY




LIGHTER DESIGN.
FREER TOMORROW.

YATLL INDUSTRY CO.,LTD



CARBON FIBER MOBILITY SCOOTER

ONLY 14KG



PREMIUM CARBON FIBER

Lightweight Strength, Premium Quality

INTEGRATED CARBON FIBER TILLER

Lightweight,
Integrated Control




ADJUSTABLE ARMREST WIDTH

Flexible Adjustment,
Personalized Comfort



INTEGRATED SUPPORT FRAME

Strong Support,
Lightweight Design



STABLE FOUR- WHEEL DESIGN

Smooth,
Stable and
Comfortable



QUICK FOLDING & DUAL BATTERIES

Easy Storage,
Flexible Charging

